

Dry electrodes for EEG

EEG electrodes are produced in 2 types:

- Dry active electrode MCScap-DrA1
- Dry passive electrode MCScap-DrP1

MCScap-DrA1 is polyurethan hermetic body with PCB of amplifier and thin and durable cable.

Amplifier use for converting high input impedance from dry sensor to low output impedance.

This solution dramatically decreases moving cable noise and external interference to cable.

Dry EEG electrodes are intended for use with soft dry sensors SoftPulse™ produced by Datwyler.

Sensors SoftPulse[™] can be in 3 types:

- Flat (for skin without hair)
- Medium with 5 mm long pins (more universal)
- Long with 8 mm long pins (for thick hair)

More information about sensors:

https://datwyler.com/technology-innovation/advanced-technologies/strategic-fields/wearables-digital-health/softpulse





Electrodes MCScap-DrA1 and MCScap-DrP1 can be preinstalled in MCScap DRY electrode cap with common connector.

Cap can have different quantity of electrodes and different types of connectors, according to requirements of EEG amplifier model and aims of application.

Specification

Active EEG electrode	
Amplification	1 ±0.0003%
Passband	0 150 kHz @ -3 dB
Input current	10 pA max
Input capacitance	7 pF typ
Input offset voltage	5 μV typ, 25 μV max
Voltage noise	0.8 μV p-p @ 0.1-10 Hz
Output impedance	27 Ohm 1%
Supply voltage	monopolar 2.2 5.5 V or bipolar ±1.1 ±2.75 V
Power supply rejection ratio	105 dB
Quiescent current	1 mA max
Cable outlet diameter	1.3 mm
Cable length	up to 2000 mm
Cable color	black
Body size	D15 x 8 mm
Body weight	< 0,8 g (up to 1,25 g with sensor)
Ingress protection	IP67 with sensor installed
Passive EEG electrode	
Cable outlet diameter	1.5 mm
Cable length	up to 2000 mm
Cable color	black
Body size	D15 x 8 mm
Body weight	< 0,8 g (up to 1,25 g with sensor)
Ingress protection	IP67 with sensor installed