



## UNCOMPROMISING EEG CAP FOR LONG-TERM MONITORING DEC NET COMBINES QUALITY OF RECORDING WITH PATIENT COMFORT UP TO 72 HOURS

Disposable EEG cap with 23 AG/AgCl electrodes according international 10-20 system.

Electrode positions: Fp1, Fp2, F3, F4, F7, F8, Fz, Cz, C3, C4, Pz, T3, T4, P3, P4, A1, A2, T5, T6, O1, O2, GND, REF.

Connection of the cap to EEG amplifier via reusable adapter with DB25 or touchproof connectors.

## Main advantages:

- √ 100% solution to the problem of cross infection
- ✓ Easy and quick mounting of the cap
- ✓ Light, breath and comfortable cap for patient

## Disposable cap allows to save much time:

- no need place of each electrodes
- √ no need preparation & special skills
- √ no need clean & dry
- no need disinfection

## Sizes:

L 54 - 60 cm 21" - 23 5/8"

**M** 48 - 54 cm 19" - 21"

**S** 42 - 48 cm 16 1/2" - 19"

Medical Computer Systems Ltd. Passage 4922, 4-2, Zelenograd, Moscow, 124460, Russia Phone: +7 495 913 31 94 E-mail: mks@mks.ru Internet: www.mks.ru Internet shop: www.mcscap.com

