



Size	Head circumference	Colored cap Size is identified by textile color	Black cap Size is identified by seam color	Chin / chest belt size	Recommended for *
XL	60-66 cm	Green 	Green 	Adult	Adults
XL/L	57-63 cm	Green + blue 	Green + blue 		Adults
L	54-60 cm	Blue 	Blue 		Adults (most)
L/M	51-57 cm	Blue + red 	Blue + red 		Teenagers, adults
M	48-54 cm	Red 	Red 	Children	Children, teenagers
M/S	45-51 cm	Red + yellow 	Red + yellow 		Children up to 5 years
S	42-48 cm	Yellow 	Yellow 		Children up to 2 years
S/XS	39-45 cm	Yellow + green 	Yellow + green 		Children up to 1 year
XS	36-42 cm	Green 	Green 	Infants	Babies up to 4 months
Inf I	32-36 cm	Blue 	Blue 		Infants up to 1st month
Inf II	28-32 cm	Red 	Red 		Newborns, premature baby
Inf III	24-28 cm	Yellow 	Yellow 		Newborns, premature baby

* Recommendations are based on WHO anthropological data <http://www.who.int/childgrowth> and ISO/TR 7250- 2:2010. This information is a reference. Before examination, measure the patient's head circumference and select a cap of the appropriate size. Choosing the correct cap size is important for correct EEG recording and patient comfort.