

# How to use 10-20 tape?

**1.**

Mark all positions by pencil

Fix at **Nasion** point

Fix at **Inion** point

head size 48-60 cm 19"-24"

MCSmap accessories for EEG

**2.**

Fix at **Right Pre-Auricular** point

Fix at **Left Pre-Auricular** point

MCSmap accessories for EEG

**3.**

Fix at **Fpz** point

Fix at **Oz** point

**4.**

Fix at **Fpz** point

Fix at **Oz** point

MCSmap accessories for EEG