



EP/ERP package

Evoked Potentials / Event-Related Potential



P300
ERP component elicited in the process of decision making



VN400
ERP component elicited to visual inclusion of meaningful words in sentence contexts



GONOGO
test of variation of decisions about spatial complexity or functional connectivity of visual objects by EEG




Mathematics
test of brain's reaction in the process of solving mathematical problems



ECPT
test to assess the ability to focus on constant activity



VCPT
test to assess the ability to focus on visual stimulus



VCPT test




P300 protocol



GONOGO test

Assembly diagram



Inside this box

- ASTIM acoustic stimulator with 5m USB cable
- ER headphones
- ER keyboard
- VGASENS optical trigger sensor
- trigger cable for connection to EEG amplifier NVX
- user manual
- software license

EP/ERP package is an option for NEUROVISOR clinical EEG system. EEG amplifier NVX36/NVX52 is not included in the EP/ERP package

For installation required PC with Windows 10 64 bit or higher

Original EEG accessories MCScap

For information visit internet shop www.mcscap.com
European Authorized Representative GVB geliMED
www.gvb-gelimed.de

EP/ERP package
unboxing
 **YouTube**

